

# Package: tdtr (via r-universe)

May 19, 2026

**Title** Read and Export Tucker-Davis Technologies Data from R

**Version** 0.0.5

**Description** Provides an R interface for reading Tucker-Davis Technologies (TDT) tank/block data through the Python 'tdt' package with reticulate, including Python-backed wrappers and explicit helpers for collecting data into ordinary R objects for downstream analysis.

**License** MIT + file LICENSE

**URL** <https://github.com/nimh-dsst/tdtr>,  
<https://nimh-dsst.r-universe.dev/tdtr>

**BugReports** <https://github.com/nimh-dsst/tdtr/issues>

**Encoding** UTF-8

**Roxygen** list(markdown = TRUE)

**Depends** R (>= 4.1.0)

**Imports** jsonlite, reticulate (>= 1.41.0), rlang, tibble

**Suggests** knitr, rmarkdown, testthat (>= 3.0.0), withr

**VignetteBuilder** knitr

**Config/testthat/edition** 3

**Config/roxygen2/version** 8.0.0

**Config/pak/sysreqs** libpng-dev python3

**Repository** <https://nimh-dsst.r-universe.dev>

**Date/Publication** 2026-05-19 16:15:16 UTC

**RemoteUrl** <https://github.com/nimh-dsst/tdtr>

**RemoteRef** HEAD

**RemoteSha** 6827dbcb8da4e6adff4b91025849b1cc75efc16b

## Contents

as_ranges . . . . .	2
as_tdt_block . . . . .	3
as_tibble_epocs . . . . .	3
as_tibble_stream . . . . .	4
block_info . . . . .	4
collect_block . . . . .	5
collect_epocs . . . . .	5
collect_stream . . . . .	6
epoc . . . . .	7
epoc_filter_py . . . . .	7
epoc_names . . . . .	8
epocs . . . . .	8
is_tdt_block . . . . .	9
new_tdt_block . . . . .	9
profile_tdt_memory . . . . .	10
ranges_from_epocs . . . . .	11
read_block . . . . .	11
read_block_py . . . . .	12
read_epocs_csv . . . . .	13
read_sev_py . . . . .	14
read_stream_csv . . . . .	15
stream . . . . .	15
stream_names . . . . .	16
streams . . . . .	16
tdt_available . . . . .	17
tdt_config . . . . .	17
tdtr_example_block_path . . . . .	18
validate_tdt_block . . . . .	18
<b>Index</b>	<b>19</b>

---

as_ranges	<i>Convert time windows to Python tdt ranges</i>
-----------	--

---

### Description

Python tdt expects a 2 x N numeric matrix where row 1 is start time and row 2 is stop time.

### Usage

```
as_ranges(ranges)
```

### Arguments

ranges	A numeric length-2 vector, a 2 x N matrix, an N x 2 matrix, or a data frame with start and stop columns.
--------	--

**Value**

A numeric 2 x N matrix.

---

as_tdt_block	<i>Coerce to a materialized TDT block</i>
--------------	---

---

**Description**

Coerce to a materialized TDT block

**Usage**

```
as_tdt_block(x, ...)
```

**Arguments**

x	A tdt_block, tdt_block_py, or compatible object.
...	Passed to collect_block().

**Value**

A materialized tdt\_block.

---

as_tibble_epocs	<i>Return epocs/events as a tibble</i>
-----------------	--

---

**Description**

Return epocs/events as a tibble

**Usage**

```
as_tibble_epocs(x, store = NULL)
```

**Arguments**

x	A tdt_block, tdt_block_py, or compatible object.
store	Optional epoc store name.

**Value**

A tibble with store, onset, offset, and value columns.

---

as_tibble_stream	<i>Return a bounded stream table</i>
------------------	--------------------------------------

---

### Description

Return a bounded stream table

### Usage

```
as_tibble_stream(x, stream, window = NULL, downsample = NULL, max_rows = NULL)
```

### Arguments

x	A tdt_block, tdt_block_py, or compatible object.
stream	Stream name.
window	Optional length-2 time window in seconds.
downsample	Optional integer stride.
max_rows	Optional maximum rows to return.

### Value

A tibble with time, channel, and value columns.

---

block_info	<i>Return block metadata</i>
------------	------------------------------

---

### Description

Return block metadata

### Usage

```
block_info(x)
```

### Arguments

x	A tdt_block, tdt_block_py, or compatible object.
---	--

### Value

A named list.

---

collect_block	<i>Collect a Python-backed TDT block into R</i>
---------------	---

---

**Description**

Collect a Python-backed TDT block into R

**Usage**

```
collect_block(
  x,
  streams = TRUE,
  epocs = TRUE,
  snips = FALSE,
  scalars = FALSE,
  stores = NULL,
  max_bytes_warn = 500 * 1024^2,
  quiet = FALSE
)
```

**Arguments**

x	A tdt_block, tdt_block_py, or compatible object.
streams	If TRUE, collect stream data.
epocs	If TRUE, collect epoc/event data.
snips, scalars	Reserved for future collection of Python tdt stores.
stores	Optional stream names to collect.
max_bytes_warn	Warn before copying a Python array larger than this many bytes.
quiet	Suppress size warnings.

**Value**

A materialized tdt\_block.

---

collect_epocs	<i>Collect epocs/events into R</i>
---------------	------------------------------------

---

**Description**

Collect epocs/events into R

**Usage**

```
collect_epocs(x, store = NULL, as = c("tibble", "list"))
```

**Arguments**

x	A tdt_block, tdt_block_py, or compatible object.
store	Optional epoc store name. If NULL, all epocs are collected.
as	Return shape: "tibble" or "list".

**Value**

A tibble or named list of tibbles.

---

collect_stream	<i>Collect one stream into R</i>
----------------	----------------------------------

---

**Description**

Collect one stream into R

**Usage**

```
collect_stream(
  x,
  store,
  as = c("matrix", "numeric", "list"),
  include_time = FALSE,
  max_bytes_warn = 500 * 1024^2,
  quiet = FALSE
)
```

**Arguments**

x	A tdt_block, tdt_block_py, tdt_sev_py, or compatible object.
store	Stream store name. The sanitized Python key or original TDT store name may be used.
as	Return shape: "matrix", "numeric", or "list".
include_time	If TRUE, return a list with data, time, and stream metadata.
max_bytes_warn	Warn before copying a Python array larger than this many bytes.
quiet	Suppress size warnings.

**Value**

A matrix, numeric vector, or stream metadata list.

---

epoc *Return one epoc/event store*

---

### Description

Return one epoc/event store

### Usage

```
epoc(x, store)
```

### Arguments

x	A tdt_block, tdt_block_py, or compatible object.
store	Epoc store name.

### Value

One epoc object or tibble subset.

---

epoc\_filter\_py *Filter TDT epocs through Python tdt*

---

### Description

Filter TDT epocs through Python tdt

### Usage

```
epoc_filter_py(
  data,
  epoc,
  values = NULL,
  modifiers = NULL,
  t = NULL,
  tref = FALSE,
  keepdata = TRUE,
  ...
)
```

### Arguments

data	A tdt_block_py wrapper or raw Python block object.
epoc	Epoc store name.
values, modifiers, t, tref, keepdata	Arguments passed to Python tdt.epoc_filter().
...	Additional keyword arguments passed to Python tdt.epoc_filter().

**Value**

A `tdt_block_py` wrapper.

---

<code>epoc_names</code>	<i>Return epoc/event store names</i>
-------------------------	--------------------------------------

---

**Description**

Return epoc/event store names

**Usage**

```
epoc_names(x)
```

**Arguments**

`x` A `tdt_block`, `tdt_block_py`, epoc container, or compatible object.

**Value**

A character vector of epoc store names.

---

<code>epocs</code>	<i>Return the epoc/event container or rows</i>
--------------------	--

---

**Description**

Return the epoc/event container or rows

**Usage**

```
epocs(x, store = NULL)
```

**Arguments**

`x` A `tdt_block`, `tdt_block_py`, or compatible object.

`store` Optional epoc store name.

**Value**

For materialized blocks, a tibble. For Python-backed blocks with no store, the live Python epocs container.

---

is_tdt_block	<i>Test whether an object is a materialized TDT block</i>
--------------	---

---

**Description**

Test whether an object is a materialized TDT block

**Usage**

```
is_tdt_block(x)
```

**Arguments**

x	Object to test.
---	-----------------

**Value**

TRUE if x inherits from tdt\_block, otherwise FALSE.

---

new_tdt_block	<i>Create a materialized TDT block</i>
---------------	--

---

**Description**

Create a materialized TDT block

**Usage**

```
new_tdt_block(info = list(), streams = list(), epocs = NULL)
```

**Arguments**

info	Named list of block metadata.
streams	Named list of stream objects.
epocs	Epoc/event table. If NULL, an empty tibble is used.

**Value**

A tdt\_block object.

---

profile\_tdt\_memory      *Profile memory movement for a TDT read workflow*

---

## Description

profile\_tdt\_memory() runs a small read workflow and records elapsed time, R allocations reported by [Rprofmem\(\)](#), and Python allocations reported by Python tracemalloc when available. Use reader arguments such as store, channel, t1, t2, evtype, and ranges to test a bounded workflow before scaling up to a larger block.

## Usage

```
profile_tdt_memory(
  block_path,
  ...,
  stream = NULL,
  events = TRUE,
  summarize = TRUE,
  quiet = TRUE
)
```

## Arguments

block_path	Path to a TDT block directory.
...	Arguments passed to <a href="#">read_block()</a> .
stream	Optional stream store to collect after the read. Leave NULL to profile metadata and events without copying stream data into R.
events	If TRUE, collect epochs/events after the read.
summarize	If TRUE, call <a href="#">summary()</a> on the Python-backed block.
quiet	Suppress collection size warnings.

## Details

The numbers are diagnostics, not portable benchmarks. They are most useful for comparing alternative reads on the same machine and Python environment.

## Value

A tibble with one row per profiled step.

---

ranges_from_epocs	<i>Build time ranges from epoc/event onsets</i>
-------------------	---

---

**Description**

Build time ranges from epoc/event onsets

**Usage**

```
ranges_from_epocs(epocs, pre, post, onset_col = "onset", drop_negative = TRUE)
```

**Arguments**

epocs	A data frame with an onset column.
pre	Start offset in seconds relative to onset. Use a negative value for time before onset.
post	Stop offset in seconds relative to onset.
onset_col	Name of the onset column.
drop_negative	If TRUE, drop ranges with negative start times.

**Value**

A numeric 2 x N matrix.

---

read_block	<i>Read a TDT block through Python tdt</i>
------------	--

---

**Description**

Read a TDT block through Python tdt

**Usage**

```
read_block(block_path, ..., collect = FALSE)
```

**Arguments**

block_path	Path to a TDT block/tank directory.
...	Extra keyword arguments passed to <code>tdt.read_block()</code> .
collect	If FALSE, return a Python-backed <code>tdt_block_py</code> wrapper. If TRUE, collect into a materialized R <code>tdt_block</code> .

**Details**

read\_block() normalizes read-time store filters against the block header, so the stream names returned by stream\_names() can also be used for filtering. This matters for TDT stores whose original names are not valid Python identifiers. For example, Python tdt filters a store originally named 465A with store = "465A", but the returned stream is exposed under the sanitized name \_465A. read\_block() accepts either spelling. Use read\_block\_py() when you need Python tdt.read\_block() store matching exactly as implemented upstream.

**Value**

A tdt\_block\_py wrapper or a materialized tdt\_block.

---

read_block_py	<i>Read a TDT block through Python tdt</i>
---------------	--

---

**Description**

This is the explicit Python-backed compatibility wrapper for tdt.read\_block(). It keeps Python objects live and does not copy stream arrays into R unless a collection helper is called later.

**Usage**

```
read_block_py(
  block_path,
  bitwise = "",
  channel = 0,
  combine = NULL,
  headers = 0,
  nodata = FALSE,
  ranges = NULL,
  store = "",
  t1 = 0,
  t2 = 0,
  evtype = NULL,
  verbose = 0,
  sortname = "TankSort",
  export = NULL,
  scale = 1,
  dtype = NULL,
  outdir = NULL,
  prefix = NULL,
  outfile = NULL,
  dmy = FALSE,
  noepocauto = FALSE,
  ...
)
```

**Arguments**

block\_path      Path to a TDT block/tank directory.  
 bitwise, channel, combine, headers, nodata, ranges, store, t1, t2, evtype,  
 verbose, sortname, export, scale, dtype, outdir, prefix, outfile, dmy,  
 noepocauto      Arguments passed to Python `tdt.read_block()`.  
 ...              Additional keyword arguments passed to Python `tdt.read_block()`.

**Details**

`read_block_py()` passes `store` directly to Python `tdt.read_block()`. For stores whose returned names are sanitized by Python, such as `_465A`, pass the original TDT store ID (`465A`) or use [read\\_block\(\)](#) for `tdtr`'s header-based normalization.

**Value**

A `tdt_block_py` wrapper.

---

<code>read_epocs_csv</code>	<i>Read an epoc/event CSV export</i>
-----------------------------	--------------------------------------

---

**Description**

Read an epoc/event CSV export

**Usage**

```
read_epocs_csv(path, ...)
```

**Arguments**

path            CSV file path.  
 ...            Passed to [utils::read.csv\(\)](#).

**Value**

A tibble with `store`, `onset`, `offset`, and `value` columns.

---

`read_sev_py`*Read TDT SEV data through Python tdt*

---

## Description

Read TDT SEV data through Python tdt

## Usage

```
read_sev_py(  
    sev_dir,  
    channel = 0,  
    event_name = "",  
    t1 = 0,  
    t2 = 0,  
    fs = 0,  
    ranges = NULL,  
    verbose = 0,  
    just_names = FALSE,  
    export = NULL,  
    scale = 1,  
    dtype = NULL,  
    outdir = NULL,  
    prefix = NULL,  
    ...  
)
```

## Arguments

`sev_dir` Path to a SEV directory or file.  
`channel`, `event_name`, `t1`, `t2`, `fs`, `ranges`, `verbose`, `just_names`, `export`, `scale`,  
`dtype`, `outdir`, `prefix` Arguments passed to Python `tdt.read_sev()`.  
... Additional keyword arguments passed to Python `tdt.read_sev()`.

## Value

A `tdt_sev_py` wrapper.

---

read_stream_csv	<i>Read a stream CSV export</i>
-----------------	---------------------------------

---

**Description**

Read a stream CSV export

**Usage**

```
read_stream_csv(path, fs, name, channels = NULL, t0 = 0, ...)
```

**Arguments**

path	CSV file path.
fs	Sampling frequency.
name	Stream name.
channels	Optional channel labels.
t0	Stream start time in seconds.
...	Passed to <code>utils::read.csv()</code> .

**Value**

A materialized stream object.

---

stream	<i>Return one stream</i>
--------	--------------------------

---

**Description**

Return one stream

**Usage**

```
stream(x, name)
```

**Arguments**

x	A <code>tdt_block</code> , <code>tdt_block_py</code> , <code>tdt_sev_py</code> , or compatible object.
name	Stream name. The sanitized Python key or original TDT store name may be used.

**Value**

A stream object. For Python-backed blocks this is a live Python stream object and does not copy stream data into R.

---

stream_names	<i>Return stream names</i>
--------------	----------------------------

---

**Description**

Return stream names

**Usage**

```
stream_names(x)
```

**Arguments**

x            A tdt\_block, tdt\_block\_py, tdt\_sev\_py, stream container, or compatible object.

**Value**

A character vector of stream names.

---

streams	<i>Return the stream container</i>
---------	------------------------------------

---

**Description**

Return the stream container

**Usage**

```
streams(x)
```

**Arguments**

x            A tdt\_block, tdt\_block\_py, tdt\_sev\_py, or compatible object.

**Value**

A stream container. For Python-backed objects this may be a live Python object and does not copy stream arrays into R.

---

tdt_available	<i>Check whether Python tdt is available</i>
---------------	--

---

**Description**

Check whether Python tdt is available

**Usage**

```
tdt_available(initialize = TRUE)
```

**Arguments**

`initialize`      If TRUE, allow reticulate to initialize Python while checking for tdt.

**Value**

TRUE if reticulate can find Python tdt, otherwise FALSE.

---

tdt_config	<i>Show Python configuration for tdt</i>
------------	--

---

**Description**

Show Python configuration for tdt

**Usage**

```
tdt_config(initialize = FALSE)
```

**Arguments**

`initialize`      If TRUE, include `reticulate::py_config()` details. This may initialize Python.

**Value**

A list with reticulate, Python, and Python tdt availability fields.

---

tdtr\_example\_block\_path

*Locate the packaged TDT example block*

---

**Description**

tdtr includes a small raw TDT block from the official TDT example data so examples, tests, and vignettes can exercise the reticulate-backed reader. The files are stored in the same block layout used by TDT, under inst/extdata.

**Usage**

```
tdtr_example_block_path(mustWork = TRUE)
```

**Arguments**

mustWork            If TRUE, error when the example block is not available.

**Value**

A single path to the packaged Subject1-211115-094936 block.

---

validate\_tdt\_block

*Validate a materialized TDT block*

---

**Description**

Validate a materialized TDT block

**Usage**

```
validate_tdt_block(x)
```

**Arguments**

x                    Object to validate.

**Value**

x, invisibly.

# Index

as\_ranges, 2  
as\_tdt\_block, 3  
as\_tibble\_epocs, 3  
as\_tibble\_stream, 4

block\_info, 4

collect\_block, 5  
collect\_epocs, 5  
collect\_stream, 6

epoc, 7  
epoc\_filter\_py, 7  
epoc\_names, 8  
epocs, 8

is\_tdt\_block, 9

new\_tdt\_block, 9

profile\_tdt\_memory, 10

ranges\_from\_epocs, 11  
read\_block, 11  
read\_block(), 10, 13  
read\_block\_py, 12  
read\_block\_py(), 12  
read\_epocs\_csv, 13  
read\_sev\_py, 14  
read\_stream\_csv, 15  
Rprofmem(), 10

stream, 15  
stream\_names, 16  
stream\_names(), 12  
streams, 16  
summary(), 10

tdt\_available, 17  
tdt\_config, 17  
tdtr\_example\_block\_path, 18

utils::read.csv(), 13, 15  
validate\_tdt\_block, 18